


Servo-u

– short guide

Setting up


1. Plug in the power cord.
2. Open the hatch on the side and switch the ventilator to the on position.
NOTE: When switching on the Servo-u, you need to first pull the ON/OFF switch downwards.
3. Connect the air and oxygen hoses.
4. Lock the wheels on the cart.
5. Start the **PRE-USE CHECK**. Follow on-screen instructions. (You need the test tube during the Pre-Use Check). The patient circuit test is included in the Pre-Use Check.
6. Choose **PATIENT CATEGORY**: adult, pediatric or neonatal (optional).
7. Choose **VENTILATION TYPE**: invasive or non-invasive (optional).
8. Calculate Predicted Body Weight (PBW) automatically by tapping **PBW** and enter patient gender and height. (For pediatric and neonatal patient category it is the patient's weight). Tidal Volume in ml/kg will now be monitored.
9. Tap **VENTILATION MODES**. (If you tap and hold on the mode you will get more information).
10. Choose mode and adjust the settings.
11. Some settings have dynamic images.
12. The safety scale can be expanded by tapping **+**.
13. Change the alarm limits in alarm window.
14. Tap **START VENTILATION**.


Refer to User's manual for detailed information



NOTE: Alarms can be in one of three colors depending on priority:

- High priority**
- Medium priority**
- Low priority**

 Alarms audio paused – Active alarms can be silenced for two minutes by tapping.

 Pre-silence alarms – To silence most alarms for two minutes, tap when no alarms are active.



Tap on the active alarm and you can adjust alarm limits directly.



Current measured value

Autoset Alarm limits in controlled modes



Alarm volume

Stop ventilation: tap and hold **STANDBY**

Disconnection

Maneuvers:

Open Lung tool (option), Nebulization, Static Measurements, Manual Breath, O₂ boost level and Pes Catheter Positioning (option)

Different views



CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people, today and tomorrow.



Manufacturer - Maquet Critical Care AB - Röntgenvägen 2 - SE-171 54 Solna - Sweden - +46 (0)10 335 73 00
Sales Office - Getinge - 1 Geoffrey Way - Wayne, NJ 07470 - USA

www.getinge.com