




T-DOC v16 Utilizing User Fields

1. User Fields can be applied to and edited in most data set types in T-DOC Admin. For this guide we will utilize User Fields in Products. Click **Products**.


Note: If you do not have Products on your fast Menu refer to **Customize Fast Menu** T-DOC Quick Guide.




2. Click the **user fields** tab. To edit or add new fields click .



3. To create a new user field click an unused user field number and rename.



4. Click the Type drop box to select the type of user field.



5. Enter values and click ok.

The screenshot shows a dialog box titled "List of user fields". On the left, there is a list of 20 user fields. Fields 2 through 7 are checked, and field 8, "Vendor Tray", is selected. Fields 9 through 20 are unchecked. On the right, the configuration for the selected field is shown: "Name" is "Vendor Tray", "Type" is "Checkbox", and "Checkbox values" are "No" and "Yes". There is an unchecked checkbox for "Enable undefined value". At the bottom, there are "OK" and "Cancel" buttons.

6. The new user field will appear under the user field tab.

The screenshot shows a software interface with several tabs: "General", "Packing", "Wash / ST", "Price", "Triggers", "Indicators", "Timing", "Identical", and "User field". The "User field" tab is active. It contains a list of user fields: "OR Stock Placement (2)", "Consignment (Enter Y or N) (3)", "Intensity Value 2 (4)", "Intensity Value 3 (5)", "Implant (Enter Y or N) (6)", "Consignment 2 (7)", and "Vendor Tray (8)". The "Vendor Tray (8)" field is currently unchecked. At the bottom, there are "Edit", "Save", and "Cancel" buttons.

Note: User fields can be utilized in advanced search. Refer to Advanced Search T-DOC Quick Guide.

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. • Copyright 2022 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Sales Office, USA • 1 Geoffrey Way • Wayne, NJ 07470 • USA • Tdocsupport.us@getinge

www.getinge.com

